

Daily Treatment Patient Verification Chart — A0

| | | | |
|------------------|--------------------|----------------------|-------------------|
| Patient Name | Alice Brown | Patient ID | ABM22334 |
| DOB | 1960-04-12 | Radiation Oncologist | Dr. Michael Evans |
| Treatment Site | Femur (Palliation) | Plan/Course ID | ABM22334_FemurPal |
| Fraction | 1 of 5 | Dose/Fx (cGy) | 200 |
| Total Dose (cGy) | 1000 | Imaging Type | kV-kV |
| Bolus | None | | |

Therapist Alerts / Special Instructions

- kV-kV setup pair required prior to treatment.

Planned Fields / Steps

| Field/Step | Technique | Energy (MV) | MU | Gantry | Coll | Jaws (cm) |
|------------------|-----------|-------------|-----|--------|------|----------------------------------|
| AP Femur | STATIC | 6 | 155 | 0 | 0 | X1 -7.5, X2 7.5, Y1 -5.0, Y2 5.0 |
| kV-kV Setup Pair | Imaging | 0 | 0 | 0 | 0 | X1 0.0, X2 0.0, Y1 0.0, Y2 0.0 |
| PA Femur | STATIC | 6 | 160 | 180 | 0 | X1 -7.5, X2 7.5, Y1 -5.0, Y2 5.0 |

Use this chart to verify that all on-screen parameters match prior to Beam On.