

# Daily Treatment Patient Verification Chart — B0

Patient Name	George Harris	Patient ID	GHS12345
DOB	1955-08-03	Radiation Oncologist	Dr. Sarah Lewis
Treatment Site	Skeletal Spine	Plan/Course ID	GHS12345_SkeletalSpine
Fraction	1 of 10	Dose/Fx (cGy)	300
Total Dose (cGy)	3000	Imaging Type	kV-kV
Bolus	None		

### Therapist Alerts / Special Instructions

- kV-kV setup pair required prior to treatment.

### Planned Fields / Steps

Field/Step	Technique	Energy (MV)	MU	Gantry	Coll	Jaws (cm)
AP Spine	STATIC	6	180	0	0	X1 -8.0, X2 8.0, Y1 -6.0, Y2 6.0
kV-kV Setup Pair	Imaging	0	0	0	0	X1 0.0, X2 0.0, Y1 0.0, Y2 0.0
PA Spine	STATIC	6	185	180	0	X1 -8.0, X2 8.0, Y1 -6.0, Y2 6.0

Use this chart to verify that all on-screen parameters match prior to Beam On.